

Manual Treatment Tables

Versatility Series

The Hamilton Medical Versatility series of fixed-height treatment tables feature a reinforced cross-base frame and outrigger base design which maximizes patient comfort and maintains optimal efficiency of floor space.

Treatment Flat Top #6K9

All features of the Versatility Standard plus:

- Reinforced cross-base frame design
- Load: 500 lbs. (226.8kg) evenly distributed
- Outrigger design with adjustable leveling feet



Treatment Adjustable Back #6K8

All features of the Versatility Standard plus:

- Reinforced cross-base frame design
- Gas-spring assisted back
- Infinite position back
- Load: 500 lbs. evenly distributed
- Outrigger design with adjustable leveling feet



Colors

- Upholstery available in 5 standard vinyl colors
- Custom colors available

**Due to the limitations of the printing process, actual colors may vary slightly*



Specifications

- Height: 30in. (76.2cm)
- Width: 29in. (73.7cm)
- Length: 72in. (182.9cm) or 78in. (198.1cm)
- Backrest length: 33in. (83.82cm)
- Weight: 300lbs. (136kg)
- Load: 500lbs. (226.8kg)

Standard Features

Available on all Versatility models

- Powder-coated steel frame (chameleon)
- Seamless top, easy to clean, stain resistant upholstery
- Comfortable, high-density foam
- Ster-O-Sheet paper roll holder
- Heavy duty steel base for durability and easy maintenance

Options

- All-steel two-drawer unit
- Drawer Size: 19 3/8 x 17 x 4in. (49.2 x 43.2 x 10.2 cm)
- Extended 78 in. length top (78")
- Pneumatic elevating one-touch backrest

Room Accessories

- Step stool
- Step stool with hand rail
- Examination stools
- Guest/patient chair
- Treatment cabinet
- Waste management receptacle

Hamiltonmedical

Phone: 888-871-5600 | Fax: 843-871-2007

www.hamiltonmed.com | sales@hamiltonmed.com | 128 Berkeley Circle | Summerville, SC 29483