

Power Exam Tables

Medallion Series

The Hamilton Medical Medallion series of power exam tables and accessories are designed to be our economical offering for the basic needs of the power healthcare examination suite.



Medallion II
with premium top
and drawer options



Medallion I
with drawer unit



Medallion I
with Bierhoff crutches



Medallion
ExamCenter

Medallion I #3K5

All of the Medallion Standard features plus:

- Gas-spring assisted tilt backrest

Medallion II #3K52

All of the Medallion Standard features plus:

- Powered back position with quiet electro-mechanical drive

Medallion ExamCenter

All features of the Medallion I or the Medallion II plus:

- 26in. top
- Accessory rail

Colors

- Upholstery available in 5 standard vinyl colors
- Custom colors available

**Due to the limitations of the printing process, actual colors may vary slightly*



Specifications

- Height: 23 - 37in. (58.4 - 94cm)
- Width: 29in. (73.7cm)
- Length: 72in. (182.9cm)
- Weight: 415lbs. (188kg)
- Load: 400lbs. (181.4kg)

Standard Features

Available on all Medallion models

- Powder-coated steel frame (chameleon)
- Integrated treatment slide-basin leg rest system
- Ster-O-Sheet paper roll holder
- Seamless top, easy to clean, stain resistant upholstery
- Pneumatic elevating one-touch backrest
- Powered height position with quiet electromechanical drive
- Foot control and duplex electrical outlet
- Pelvic tilt
- Stirrups

Options

- Premium top
- Instrument warmer
- Bierhoff crutches
- Stainless steel treatment/drain pans
- All steel two-drawer unit
- Drawer Size: 19 3/8 x 17 x 4in. (49.2 x 43.2 x 10.2 cm)

Room Accessories

- Examination stools
- Guest/patient chairs
- Treatment cabinet
- Waste management receptacle

ExamCenter Accessories

Configuration:

- 26in. width top
- Accessory rail

Plus Your Choice of Accessories:

- Drape shield
- IV pole
- Padded armboard
- Armrests
- Restraint straps

Hamiltomedical

Phone: 888-871-5600 | Fax: 843-871-2007

www.hamiltonmed.com | sales@hamiltonmed.com | 128 Berkeley Circle | Summerville, SC 29483