

Power Treatment Tables

Versatility Series

The Hamilton Medical Versatility series of variable-height treatment tables feature a table base design which maximizes patient comfort and maintains optimal efficiency of floor space.

Hi-Lo Flat Top

All features of the Versatility Standard plus:

- Pedestal base with quiet electromechanical drive
- Foot control and duplex electrical outlet
- Load: 400 lbs. evenly distributed



Hi-Lo
Power Flat Top

Hi-Lo Adjustable Back

All features of the Versatility Hi-Lo Flat Top plus:

- Powered height position with quiet electro-mechanical drive
- Infinite-position back setting
- Pneumatic infinite position adjustable backrest



Hi-Lo
Power Tilt Top

Colors

- Upholstery available in 5 standard vinyl colors
- Custom colors available

**Due to the limitations of the printing process, actual colors may vary slightly*



Specifications

- Height: 30in. (76.2cm)
- Width: 29in. (73.7cm)
- Length: 72in. (182.9cm) or 78in (198.1cm)
- Weight: 415lbs. (188kg)
- Load: 400lbs. (181.4kg)

Standard Features

Available on all Versatility models

- Seamless top, easy to clean, stain resistant upholstery
- Comfortable, high-density foam
- Ster-O-Sheet paper roll holder
- Powder-coated steel frame (chameleon)
- Heavy duty steel base for durability and easy maintenance

Options

- Pneumatic infinite position adjustable backrest
- One or two all-steel two-drawer units
- Drawer Size: 19 ⁵/₈ x 17 x 4in. (49.2 x 43.2 x 10.2 cm)
- Extended 78 in. length top (198.1cm)

Room Accessories

- Step stool
- Step stool with hand rail
- Examination stools
- Guest/patient chair
- Treatment cabinet
- Waste management receptacle

Hamiltomedical

Phone: 888-871-5600 | Fax: 843-871-2007

www.hamiltonmed.com | sales@hamiltonmed.com | 128 Berkeley Circle | Summerville, SC 29483