

Seating

Blood Drawing Chairs

Heavy-duty steel base with two, height adjustable chair armrests and durable, reinforced seat for patient comfort and security. Easy to clean surface, two-position armrest for patient immobility, easy access and superior safety.

Traditional Phlebotomy Chair #7K12

Specifications:

- Seat Height: 20.5in. (52cm)
- Seat Width: 18in. (46cm)
- Overall Height: 36in. (91.4cm)
- Overall Width: 29in. (73.7cm)
- Armrest Adj. Height: 28in.-38in. (71.1-96.5cm)
- Overall Depth: 20.5in. (52cm)
- Weight: 30 lbs (14kg)
- Weight Capacity: 300lbs (136kg)

Features:

- Heavy Duty Base
- Two Height Adjustable Armrests
- Durable Fiberglass-Reinforced One Piece Black Seat/Back Access
- Easy-To-Clean



Padded Phlebotomy Chair #7K73

Specifications:

- Seat Height: 22.5in. (57.2cm)
- Seat Width: 18.75in. (48cm)
- Overall Height: 44.5in. (113cm)
- Overall Width: 29in. (73.7cm)
- Armrest Adj. Height: 28.5in.-38.5in. (72.4-97.8cm)
- Overall Depth: 27in. (69cm)
- Weight: 48 lbs (21.8kg)
- Weight Capacity: 300 lbs (136kg)

Features:

- Two Infinitely Adjustable Armrests
- Padded Armrests
- Adjustable Fold-Up Padded Armrest for Safety and Easy Access
- Padded Seat and Back
- Non-Marring Feet Levelers



Bariatric Phlebotomy Chair #7K75

Specifications:

- Seat Height: 22.5in. (57.2cm)
- Seat Width: 24.25in. (62cm)
- Overall Height: 36.5in. (92.7cm)
- Overall Width: 34in. (86.4cm)
- Armrest Adj. Height: 28.5in.-38.5in. (72.4-97.8cm)
- Overall Depth: 25in. (64cm)
- Weight: 53 lbs (24kg)
- Weight Capacity: 450 lbs (204kg)

Features:

- Two Infinitely Adjustable Armrests
- Padded Armrests
- Adjustable Fold-Up Padded Armrest for Safety and Easy Access
- Extra Large Padded Seat and Back
- Non-Marring Feet Levelers



Colors

- Upholstery available in 5 standard vinyl colors
- Custom colors available

**Due to the limitations of the printing process, actual colors may vary slightly*



Hamiltomedical

Phone: 888-871-5600 | Fax: 843-871-2007

www.hamiltonmed.com | sales@hamiltonmed.com | 128 Berkeley Circle | Summerville, SC 29483