

# Manual Treatment Tables

## V2 Series

The Hamilton Medical Basic series of fixed-height flat and adjustable back treatment tables feature a sturdy base design which maximizes patient comfort while maintaining optimal efficiency of floor space.

### V2 Standard #6K50F

### V2 Standard Shelf #6K50FS

### V2 Standard Drawer #6K50FD



V2 Flat Top  
#6K50F



V2 Flat Top  
with optional shelf unit  
#6K50FS



V2 Flat Top  
with optional drawer unit  
#6K50FD

### V2 Adjustable Back #6K57A

### V2 Adjustable Back with Shelf #6K57AS

### V2 Adjustable Back with Drawer #6K57AD



V2 Tilt Top  
#6K57A



V2 Tilt Top  
with optional shelf unit  
#6K57AS



V2 Tilt Top  
with optional drawer unit  
#6K57AD

## Colors

- Upholstery available in 5 standard vinyl colors
- Custom colors available

*\*Due to the limitations of the printing process, actual colors may vary slightly*



## Specifications

- Height: 30in. (76.2cm)
- Width: 28in. (71.1cm)
- Length: 72in. (182.9cm)
- Weight: 74 - 198lbs. (33.6 - 89.8kg)
- Load: 400lbs. (181.4kg)

## Standard Features

Available on all V2 models

- Seamless top, easy to clean, stain resistant upholstery
- Comfortable, medium-density foam
- Ster-O-Sheet paper roll holder
- Reinforced H-brace for superior strength and stability
- Powder-coated steel
- Adjustable leveling feet

## Options

- All-steel two-drawer units
- Full-length shelf beneath table
- Pneumatic elevating backrest
- 25in. (63.5cm) and 36in. (91.4cm) fixed-height models available

## Accessories

- Step stool
- Step stool with hand rail
- Examination stools
- Guest/patient chair
- Treatment cabinet
- Waste management receptacle

# Hamiltonmedical

Phone: 888-871-5600 | Fax: 843-871-2007

www.hamiltonmed.com | sales@hamiltonmed.com | 128 Berkeley Circle | Summerville, SC 29483